


WWW.IJOY.COM

 INTERACTIVE HEALTH

© Interactive Health[®], 2006
Patents Pending and intellectual property rights enforced.
No medical claims are warranted or implied by the use of this product.
iJoy is a registered trademark of Interactive Health.
www.interhealth.com

iJoy-EQ400-A0



iJoy | **Equalizer**
LIVE. PLAY. RELAX. ⁴⁰⁰
USE & CARE MANUAL

iJoy
LIVE. PLAY. RELAX.



Get MORE iJoy® in your life!

iJoy® | Equalizer® 400

AN EXCELLENT
CHOICE!

Congratulations on choosing the iJoy® Equalizer®—the ultimate in robotic foot massage technology.

The iJoy® Equalizer® provides a deep percussion massage for tired and aching feet. Its unique action soothes your soles, relieving stress and tension. In just a few minutes, your feet will feel relaxed and rejuvenated.

Before using the iJoy® Equalizer®, read this manual. Then sit back, relax, and iJoy!

KEEP IT

SAFE



DANGER: To reduce the risk of electrical shock: Always unplug this product from the electrical outlet immediately after using and before cleaning. Do not use while bathing or in a shower. Do not place or store the product where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid. If the product falls into water, unplug it immediately. Do not reach into the water.



WARNING:

To reduce the risk of burns, fire, electric shock or injury to persons:

Keep away from water.

Do not use this product near water contained in bathtubs, basins or other vessels.

An appliance should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.

Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.

Close supervision is necessary when this appliance is used by, on, or near children, invalids, or disabled persons.

Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.

Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to a service center for examination and repair.

Do not carry this product by the power cord or use the cord as a handle.

Always keep the power cord away from heated surfaces.

Never drop or insert any object into any opening.

Do not use outdoors.

Do not operate where aerosol (spray) products are being used or where oxygen is being administered. To disconnect, turn all controls to the OFF position, then remove plug from outlet.

Never operate on a soft surface such as a bed or couch where the air openings may be blocked.

Do not use this product when pain is present or if you are pregnant. Consult a physician with any questions.


Do not use this product when standing on a damp floor, or when any part of the body is in contact with plumbing or any similar ground.

If you have any medical problems such as Phlebitis, Diabetes, or any vascular problems in your calves/ legs/ feet, consult a physician before using this product.

GETTING STARTED



1. Plug the Equalizer® into a polarized outlet.
2. Sitting in your favorite chair, remove your shoes, then place your feet on the massage pods.




To reduce the risk of electric shock, the iJoy® Equalizer® has a polarized plug. The plug will only fit into a polarized outlet. If you cannot plug your Equalizer® in, do not attempt to modify the plug or outlet yourself. Contact an electrician.






1. Press the ON/OFF side of the mode select button to turn the Equalizer® on.
2. Press the PROGRAM side of the mode select button to toggle through the four speeds and Program Mode. (Program Mode cycles through all four speeds.)



Do not massage a single area for more than 15 minutes. If your skin feels uncomfortable, or appears irritated, turn off the unit immediately. **Do not attempt to massage any part of the body other than the base of your feet.**

Reflexology Foot Chart

Based on the ancient principle that reflex areas in the feet correspond to specific parts of the anatomy. Massaging these reflex points is essential for an active lifestyle and balanced body. Reflexology helps relieve stress and tension, stimulates circulation, and promotes the flow of nerve impulses.





TAKE CARE

Product Care and Storage

Maintenance and Care:

Should your Interactive Health™ product become dirty, wipe clean with a cloth dampened with a non-abrasive cleanser. Do not use thinner, benzene or any other solvents. The product is specially designed and constructed to be maintenance free. No periodic lubrication or servicing is required. Your product should only be serviced by an authorized service center.

Care should be taken not to overload the unit. Excessive feet pressure on the massage nodes that significantly slows or stops the motor will potentially cause damage to the unit. Foreign objects should not be inserted between the massage nodes and the unit housing.

Storage:

Your product has been carefully designed to give you many years of trouble-free service. Always unplug the unit when not in use. Additionally, if not in use for a long period of time, we recommend that the power cord be coiled and that the unit be placed in an environment free of dust and moisture.

Do not store near heat or open flame. Do not leave this product exposed to direct sunlight for long periods of time. This may cause fading or damage.

**Interactive Health®
Partners in Your Well-Being**

For decades, Interactive Health has been a leader in the field of massage products. We create state-of-the-art products that embrace the latest technology while applying the ancient healing art of massage. With today's fast-paced lifestyles, balancing your physical and mental needs requires the right mix of proper nutrition, exercise, and massage. Interactive Health® Massage Systems are vital components for achieving balance between physical and mental well-being.

To learn more about our company and our products, we invite you to visit us on the Internet at: www.interhealth.com or call our Customer Service Representatives at 800-355-2762.



HTT, recommended by professionals,
designed for you.



"HTT" is the preferred technology of the ACCO.
www.accoweb.org



Specifications*

| | |
|------------------------------|---------------------------------------|
| Use: | Household (not for institutional use) |
| Operating Voltage: | AC 110-120 V / 50-60 Hz |
| Watts: | 30 Watts |
| Number of Massage Functions: | 1 |
| PPM (Pulses Per Minute): | 1900-2600 |
| Product Weight: | 9.5 Lbs. |

No medical claims warranted or implied by the use of this product.

Limited One (1) Year Warranty

Please read the following warranty information for important disclosures before using your new Human Touch Technology™, iJoy® Equalizer®.

It is Interactive Health's mission to design, develop, manufacture and sell innovative, high quality massage products that create a better and healthier way of life. We appreciate your commitment to our mission and our products and are pleased to let you know that the Interactive Health®, iJoy® Equalizer® that you purchased comes with a limited one-year warranty. The warranty provides peace of mind that you are covered for mechanical or structural defects in workmanship within the coverage guidelines and limitations outlined below.

Coverage Guidelines:

The warranty begins on the original purchase date and consists of:

- **Factory Service (Labor):** Interactive Health® will supply packing materials, cover return shipping to Interactive Health®, repair parts covered by this warranty and labor for repair work and return shipping to purchaser when repairs are completed.
- **Parts Only (Parts):** Interactive Health® will supply repair parts that are covered by this warranty and return shipping to purchaser for repairs performed at the factory.
- **Structural (Metal Frame):** Interactive Health® will supply packing materials and repair parts covered by this warranty and return shipping to purchaser for repairs performed at the factory.

In the unlikely event that your iJoy® Equalizer® malfunctions or requires repair, please contact your retailer. For USA customers only, please contact Interactive Health® Customer Service at (800) 355-2762 between the hours of 8:30 AM and 5:00 PM PST or cs@interhealth.com for assistance. You can also send a fax to customer service at (800) 840-2940.

All warranty service requires approval and authorization by Interactive Health®. Products require a Return Merchandise Authorization Number (RMA) for any returns to the factory. For returns, the consumer is responsible for packing the product in the original box or in approved packaging that is available for sale from Interactive Health. Shipping and freight coverage under this warranty requires scheduling by Interactive Health and only includes "curb-side" pickup/delivery. Products received without a valid RMA number or one clearly marked on the box will be refused.

General Limitations:

- The warranty is only valid within the United States and Canada provided the product has been operated according to the instructions accompanying it and has not been misused or damaged in any way that is beyond the control of Interactive Health.
- All service covered by this warranty must be approved by Interactive Health and repairs performed by authorized technicians.
- The warranty applies only to the product; accessories or enhancements are not covered.
- Important exclusions apply to cosmetic items such as fabrics, woods, foams, padding, plastics covers, etc.
- Interactive Health is not responsible in any way for losses, damages, or inconvenience caused by equipment failure.
- For a complete listing of warranty coverage, exclusions, and limitations, visit www.interhealth.com.
- Warranties do not apply to rental, business, commercial, institutional, or other non-residential users.
- Separate warranty terms apply to used/refurbished products including those received out-of-box. Proof of purchase (original receipt) is required for all warranty repairs.

EXTENDED WARRANTIES:

Interactive Health® offers extended warranty programs to enhance ownership of the Interactive Health® product! Extended warranty purchases are available within 90 days of purchase of iJoy® product from either the retailer or direct from Interactive Health's Customer Service Department.

Please visit www.interhealth.com for the latest terms and pricing for extended warranties.

| (iJoy® Models) | Factory Service | Parts | Structural |
|--|-----------------|--------|------------|
| New Products (Original Packaging Only) | 90 days | 1 year | N/A |