

the art of style

Our trend-setting team of designers continually introduces new colors, new styles, and new fabrics, because we know it's not enough that our products provide the most effective Robotic Massage® on the market. To deserve a place in your environment, our chairs must be great looking as well as therapeutic.

It's our way of saying we understand what you want.

the science of comfort

Everyone likes to feel good. At Interactive Health, it's all we think about.

Our team of scientists, researchers, and engineers analyze the treatments used by massage therapists, chiropractors, and physicians, and turn those techniques into the sophisticated **Human Touch Technology®** mechanisms that animate our family of chairs. By building on our decades of experience as the world leader in massage technology, we consistently produce the most advanced Robotic Massage® systems ever devised.

Specifications

іЈоу™100 Product Description: iJoy™ Robotic Massage ® Chair Volt AC: 120 AMP/Watt: .58/70 Regulatory Certification: Recline Angle: 115/150 Degrees Massage Head Storage: Automatic Timer Shut-Off: 15 minutes

Total Massage Stroke Length: 20" Upholstery Material: Massage Robot:

Brushed Polyester/Nylon Mesh Human Touch Technology® Robotic Massage® Quad Rollers

Chair Reclined Dimensions: Chair Upright Dimensions: Weight Product Lbs: Master Carton Quantity: Master Carton Dimensions:

Master Carton Cube: Weight Shipping Lbs: Container Quantity (40 Foot): 176 Units Truckload Quantity (53 Foot): 252 Units Shipping Method:

Warranty:

Shipped From: Country of Origin: 46.5"L x 27"W x 24"H 41"L x 27"W x 31"H 75 lbs 1 Unit

31.5"L x 28"W x 21.25"H 11.7' 93 lbs

FOB Long Beach, CA, USA

1 Year Limited Warranty Massage Softening Pad, Use and Care Manual

1 INTERACTIVE HEALTH

13820 Stowe Drive San Diego, CA 92064 800-796-9656 •858-218-1320 FAX 858-218-1321 www.vitalityweb.com

© 2004 Interactive Health ® Patents Pending. Actual products and colors may vary. No medical claims are warranted or implied by the use of this product. iJoy is a trademark of Interactive Health.









meet your new best friend.



Inside every iJoy™ Robotic Massage® Chair are the "hands" of a trained massage professional – just waiting to provide you with a soothing back massage. But don't take our word for it. Sit down in one.

Our multi-patent **Human Touch**

Technology* massage system so closely replicates the therapeutic techniques used by back and spinal care professionals that it is the preferred technology of the American College of Chiropractic Orthopedists. It is designed to follow the natural contour of your spine, and move three-dimensionally in four different modes.



iJoy™100 does what you tell it. Each of four powerful massage modes can be programmed to give you the touch you prefer. The controller has three different programmed massage sessions of 15 minutes, which means you can kick back and enjoy a programmed massage whenever you feel like it.



Removable massag softening pad

Easy-to-reach



Child safety loo




Wired controller stores within easy rea

1-800-796-9656

www.vitalityweb.com





Ottoman 2.0™ sold separately



Rolling Kneadin
Feels like the gentle hands of a massage therapist, pressing firmly on one side of your spine, then the other. It relieves tension and loosens the muscles to prepare your improve to the spine of the spine of the spine of the spine of the spine. The muscles to prepare your improve the spine of the spi



Kneading
Feels like thumbs pressing firmly in a circular motion on both sides of the spine. This therapeutic motion goes deeper, lifting and stretching muscles to improve circulation.



Compression

(HTT™ exclusive)

Feels like the heels of a professional therapist's hands firmly pushing on muscles alongside your opine, and inward toward the spine itself. This mode aids in improving mobility

and posture.



Percussion
Feels like fists or the edge
of hands rapidly tapping on
your back muscles. This
is the deep muscle therapy
popular in sports medicine.
It flexes spinal joints, relieves
pressure, and invigorates

Functions:

Rolling: 1-Speed
Kneading: 1-Speed
Compression: 1-Speed
Percussion: 1-Speed

Range Control: Partial massage stroke 3" up & down for

kneading or compression

Position Manual Up/Down: Manually moves massage rollers to desired location for knowling, compression and persuscion

for kneading, compression and percussion

Backrest: Manual recline

Stop/Auto Park: Stops all massage functions and parks the massage

rollers in their highest position

Programmed Massage Sessions

Full: 15 minute program-Series of combined functions for your entire back Upper:15 minute program-Series of combined functions for your upper back Lower:15 minute program-Series of combined functions for your lower back