

the art of style

Our trend-setting team of designers continually introduces new colors, new styles, and new fabrics, because we know it's not enough that our products provide the most effective Robotic Massage® on the market. To deserve a place in your environment, our chairs must be great looking as well as therapeutic.

It's our way of saying we understand what you want.

the science of comfort

Everyone likes to feel good. At **Interactive Health**, it's all we think about.

Our team of scientists, researchers, and engineers analyze the treatments used by massage therapists, chiropractors, and physicians, and turn those techniques into the sophisticated **Human Touch Technology**® mechanisms that animate our family of chairs. By building on our decades of experience as the world leader in massage technology, we consistently produce the most advanced Robotic Massage® systems ever devised.

Specifications

iJoy™130 Product Description: iJoy™ Robotic Massage ® Chair Volt AC: AMP/Watt: .58/70 Regulatory Certification: Recline Angle: 110/170 Degrees Massage Head Storage: Top Automatic Timer Shut-Off: 15 minutes Total Massage Stroke Length: 20" Micro-suede/Nylon Mesh Upholstery Material: Massage Robot:

110/170 Degrees Weight Shipping L
Top Container Quantity
15 minutes Truckload Quantity
20" Shipping Method:
Micro-suede/Nylon Mesh
Human Touch Technology® Country of Origin:
Robotic Massage® Warranty:
Quad Rollers Included:

Chair Reclined Dimensions: 50"L x 29.3"W x 19.7"H Chair Upright Dimensions: 41.5"L x 29.3"W x 35.5"H Weight Product Lbs: 80 lbs Master Carton Quantity: 52"L x 30"W x 21"H Master Carton Dimensions: Master Carton Cube: 97 lbs Weight Shipping Lbs: Container Quantity (40 Foot): 135 Units Truckload Quantity (53 Foot): 180 Units Common Carrier FOB Long Beach, CA, USA P.R.C. 1 Year Limited Warranty

NINTERACTIVE HEALTH

Vitalityweb.com, INC 13820 Stowe Drive - San Diego, CA 92064 1-800-796-9656

www.vitalityweb.com

© 2004 Interactive Health ® Patents Pending. Actual products and colors may vary. No medical claims are warranted or implied by the use of this product. iJoy is a trademark of Interactive Health.



meet your new best friend.



Inside every iJoy™ Robotic Massage® Chair are the "hands" of a trained massage professional – just waiting to provide you with a soothing back massage. But don't take our word for it. Sit down in one.

Our multi-patent Human Touch

Technology® massage system so closely replicates the therapeutic techniques used by back and spinal care professionals that it is the preferred technology of the American College of Chiropractic Orthopedists. It is designed to follow the natural contour of your spine and move three-dimensionally in four different modes.



iJoy™130 does what you tell it. Each of four powerful massage modes can be programmed to give you the touch you prefer. The controller has three different programmed massage sessions of 15 minutes, which means you can kick back and enjoy a programmed massage whenever you feel like it.



softening pad







Retractable controller stores within easy reach

1-800-796-9656

www.vitalityweb.com











the science of comfort.

Rolling: Feels like the gentle hands of a massage therapist, pressing firmly on one side of your spine, then the other. It relieves tension and loosens the muscles to prepare your back for deeper massage.

Kneading: Feels like thumbs pressing firmly in a circular motion on both sides of the spine. This therapeutic motion goes deeper, lifting and stretching muscles to improve circulation.

Compression: (HTT™ exclusive) Feels like the heels of a professional therapist's hands firmly pushing on muscles alongside your spine, and inward toward the spine itself. This mode aids in improving mobility and posture.

Percussion: Feels like fists or the edge of hands rapidly tapping on your back muscles. This is the deep muscle therapy popular in sports medicine. It flexes spinal joints, relieves pressure, and invigorates your entire back.

Functions:

Rolling: 1-Speed 1-Speed Kneading: Compression: 1-Speed Percussion: 1-Speed

Range Control: Partial massage stroke 3" up & down for kneading

or compression

Position Manual Up/Down:

Manually moves massage rollers to desired location for kneading,

compression and percussion Manual recline

Backrest:

Stop/Auto Park:

Stops all massage functions and parks the massage rollers in their

highest position

Programmed Massage Sessions

Full: 15 minute program-Series of combined functions for your entire back

Upper:15 minute program-Series of combined functions for your

upper back

Lower:15 minute program-Series of combined functions for your

lower back